You are now the proud owner of The Bag!

InsulinSaver's self-cooling bag is easy to use and designed to reduce the worry of everyday storage.

Specifications:

Model	Typical duration at 35 °C/95 °F	Recharging Method	Max Temperature Protection
The Bag	11 hours +	Cooling naturally below 24 °C/75 °F or refrigerating (optional)	54 °C/129 °F

How does it work?

The Bag contains a special cooling gel that activates at temperatures above 25 °C/77 °F, changing from a solid to a liquid state as it absorbs heat. When fully "spent", the gel will feel soft, almost liquid. To recharge, allow The Bag to cool naturally in an air-conditioned room or place it in a refrigerator.

How do I know it's working?

When The Bag is charged, you'll feel the gel as firm and solid. If only some solid remains in the sides of the bag, use an **InsulinSaver Smart Temperature Sensor** to monitor the temperature conditions inside the bag.

Longevity

The Bag is built to last and does not degrade over time if used as intended. For any questions or concerns, contact us at info@insulinsaver.com

Manufacturer: InsulinSaver AB,
Gunneröd 609, 472 91 Stillingsön, Sweden.
Designed in Sweden. Made in China.

Weight 350g/0.8lb. 8x20x8 cm/3.2x8x3.2 inches.
Textile is 100% Polyester. Zipper is metal.
Element is PCM. Packaging is PE. Manual is paper.
Hand wash only.

Batch: Ba02C
ManV: 1.3.8

1. Place your insulin vials or pen inside The Bag as needed.
2. Close The Bag.
3. Feel confident that your contents are better protected from heat.

What does The Bag do?

The Bag is designed to help reduce heat exposure during warm conditions.

Tips for Optimal Use

- Place an **InsulinSaver Smart Temperature Sensor** in the small side pocket inside The Bag to monitor your medication. This setup can alert you if temperatures become too high or too low, giving you added peace of mind.
- Traveling by air? The Bag is airport-friendly! Ensure the gel is in a solid state before security checks.
- Fitting issues? If The Bag feels rigid on delivery, heat it to around 50 °C/122 °F for two hours. Then mold the sides gently before use.

Warning

Always follow your medication's storage requirements. Be aware of the recommended temperature guidelines. Never place the bag containing temperature-sensitive medication in a freezer. The Bag is not a medical device and does not diagnose, treat or prevent disease. This product is not intended for children under 12 years.

Want the manual
in your language?
Scan the QR-code!

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